

[INTENDED USE]

ADC conjugation kit is designed to conjugate antibody and payload through site-specific enzymatic reaction.

The kit is for research use only.

[PRINCIPLE OF THE Kit]






Fig. 1 The principle of DAR4 ADC preparation with AGLink®

Almost all monoclonal antibodies are glycosylated at Asn-297 of Fc domain. While the glycans have different isoforms, the types which typically dominate are G0F, G1F&G2F (>90%)^[1]. AGLink® ADC kit, which is based on the YConju™ Platform of Glyco-therapy Biotechnology, employs the N-linked glycosylation site for site-specific toxin conjugation.

The principle of conjugation is shown in Figure 1. Firstly, galactose is added to the terminal GlcNAc catalyzed by β-1,4-galactosyltransferase (GT) to form LacNAc. Then, the α-1,3-fucosyltransferase (FT) recognize LacNAc, and the toxin is transferred to the GlcNAc of LacNAc. So, the theoretical DAR value of the final conjugates is 4.

The conjugation procedure, which is enzyme catalyzed, is performed as one-pot by combining glycan modification and toxin conjugation simultaneously. The conjugation procedure, which is easy to operate, is performed under physiological conditions, and organic solvents are not introduced. The conjugates prepared through this method usually have high homogeneity and the binding activity of the conjugates to related antigens is not changed.

Overview of the protocol for constructing DAR4 ADC:



Content and Storage

BoxA

Name	Details	Quantity	Storage temperature
Enzyme A	GT	1 × 270μL	-20°C
Enzyme B	FT	1 × 120μL	-20°C
Cofactor A	MnCl2 (0.1 M)	1 × 40μL	-20°C
Cofactor B	MgCl2 (0.5 M)	1 × 40μL	-20°C
20 × Reaction buffer	20 × Reaction buffer	1 × 70 μL	-20°C
Substrate A	UDP-Gal (0.1 M)	1 × 70 μL	-20°C
Substrate B	GDP-Fuc-vc-PAB-MMAE (50 mM)	1 × 60 μL	-20°C

BoxB

Name	Details	Quantity	Storage temperature
Binding buffer	25mM Tris-HCl, 100 mM NaCl, pH 7.5	15ml	4°C
Elution buffer	0.1 M Glycine, pH 2.7	2.1ml	4°C
Neutralization buffer	1 M Tris-HCl, pH 8.0	500μl	4°C
Tris-HCl	250mM Tris-HCl buffer (pH7.5)	1.5ml	4°C
10xPBS	10xPBS	6ml	4°C
Desalting Spin column	Desalting Spin column, 2 mL	2 pieces	4°C
Ultrafiltration concentrator	Ultrafiltration concentrator (include 2 collection tubes), 0.5 mL	4 pieces	4 °C to 25 °C
Pierce™ spin column	Pierce™ spin column (containing 200 μL protein A resin)	2 pieces	4°C

Equipment required:

- Centrifuge for 1.5-2 mL&15-50 mL tubes (With refrigeration function)
- Incubator or water bath for 30 °C
- Rotary mixer
- UV/VIS Spectrophotometer (Thermo NanoDrop one)

Protocol for conjugation of 5mg of antibody:

A. Concentration and buffer exchange of the antibody

- Concentration step (Optional)

Note: Only needed when the start concentration of the antibody is lower than 11 mg/mL.

1. Add 500 μL of ddH2O to the concentrator and cap the device.
2. Centrifuge at 10000 × g for 5 min at 4°C.
3. Discard the solution in the concentrator and the flow-through.
4. Add the antibody solution to the concentrator.
5. Centrifuge at 10000 × g for 2-6 min at 4°C.

Note: If the antibody volume in the concentrator is not suitable after first concentration, centrifuge for an additional 2-6 min at 10000 × g for 4°C, until the appropriate volume is achieved.

6. Collect the solution in the concentrator.

- Buffer exchange with Desalting Spin column, 2 mL

1. Prepare 10 mL of 25mM Tris-HCl buffer (pH7.5) by adding 1 mL of 250mM Tris-HCl to 9 mL of ddH2O in a 15 mL tube. Vortex briefly to mix.
2. Break off the bottom closure of the Desalting Spin column. Loosen the lid (do not remove the lid).

3. Place the column in a collection tube (15 mL) and centrifuge at 1000 × g for 2 min at 4°C to remove the storage solution.
4. Discard the flow-through and place the column in the collection tube.
5. Add 1 mL of 25 mM Tris-HCl buffer (pH7.5) on top of the resin. Centrifuge the column at 1000 × g for 2 min at 4°C and discard the flow-through.
6. Repeat 25 mM Tris-HCl buffer (pH7.5) wash in step 5 two more times. Last spin for 3 min to remove excess liquid. Place the column in a new collection tube (15 mL).
7. Apply the antibody solution on top of the resin (100-400 μL).
8. Centrifuge at 1000 × g for 3 min at 4°C and collect the flow-through containing the antibody in 25 mL Tris-HCl buffer (pH7.5).

B. Antibody Fc domain N-glycan Modification & Toxin conjugation

1. Add appropriate concentrated antibody sample (5 mg) to a sterilized 1.5 mL tube. Note: The concentration of start antibody should be more than 11 mg/mL. If not, add by the maximum volume (470 μL).
2. Add 50 μL 20 × Reaction buffer to the antibody solution.
3. Add 50 μL UDP-Gal & 40 μL GDP-Fuc-vc-PAB-MMAE to the antibody solution.
4. Add 250 μL GT & 100 μL FT to the antibody solution.
5. Add 20 μL MnCl2 & 20 μL MgCl2 to the antibody solution.
6. Add appropriate ddH2O to the antibody solution to a whole volume of 1000 μL.
7. Mix the Reaction solution by carefully pipetting up and down.
8. Incubate the Reaction solution at incubator or water bath protected from light for 48 h, at 30 °C.

C. Removal of Enzyme and excess toxin reagent

- Equilibration

1. Place one Pierce™ spin column containing 200 μL protein A resin in 2 mL centrifuge tube.
2. Remove the cap in the bottom of the column, Centrifuge 1 minute at 1000 × g for 4°C to remove the storage solution.
3. Discard the liquid in the collection tube.
4. Add 500 μL of Binding Buffer to the column, Centrifuge 1 minute at 1000 × g for 4°C.
5. Discard the liquid in the collection tube.
6. Repeat steps 4 and 5 twice for a total of three washes.

- Binding of the antibody conjugate

7. Add 500 μL Reaction solution to the resin in the spin column. Cap top and bottom of the spin column.
8. Incubate the spin column assembly on a Rotary Mixer, and mix slightly at room temperature for 60 minutes.


- Wash

9. Remove the bottom caps and place the spin columns in 2 mL collection tubes. Loosen the top lids.
10. Centrifuge 1 minute at 1000 × g for 4°C; do not discard the collected solution.
11. Place the spin columns in a new 2 mL collection tubes.
12. Add 500 μL of Binding Buffer to the spin columns, centrifuge 1 minute at 1000 × g for 4°C.
13. Discard the liquid in the collection tube.
14. Repeat steps 12 and 13 twice for a total of three washes.
- Elution of purified, conjugated antibody
15. Cap the bottom and place the spin columns in a new 2 mL sterilized collection tubes
16. Add 500 μL Elution buffer to the resin in the spin column. Cap top of the spin column.
17. Incubate the spin column assembly on a Rotary Mixer, and mix slightly at room temperature for 10 minutes.
18. Remove the bottom caps and loosen the top lids.
19. Centrifuge the spin column assembly 1 min at 1000 × g for 4°C to elute the conjugated antibody.
20. Add 100 μL Neutralization buffer immediately to the elution to adjust the pH to neutral.
21. Repeat elution (steps 16-20) one more time.
22. Pool the eluted fractions.
23. Concentration and buffer exchange of the conjugated antibody to 1xPBS (Prepare 10 mL of 1xPSB (pH7.2-7.5) by adding 1mL of 10xPBS to 9mL of ddH2O in a 15 mL tube. Vortex briefly to mix.) according to Step A.
24. Store the conjugated antibody at -80°C for long preservation.

D. Measure the DAR value of the conjugated antibody (Herceptin as example)

•Measure the DAR value of the conjugated antibody (LC-MS based)

1. Calculation the MS of ADC



MS of DAR4 ADC= The MS of G1/G2(Herceptin)+[(180.16)(Galactose)-18.02(H₂O)]+ [(1808.8)(GDP-Fuco-vc-PAB-MMAE)-443.20(GDP)]*4

MS of DAR3 ADC= The MS of G1/G2(Herceptin)+[(180.16)(Galactose)-18.02(H₂O)]+ [(1808.8)(GDP-Fuco-vc-PAB-MMAE)-443.20(GDP)]*3

MS of DAR2 ADC= The MS of G1/G2(Herceptin)+[(180.16)(Galactose)-18.02(H₂O)]+ [(1808.8)(GDP-Fuco-vc-PAB-MMAE)-443.20(GDP)]*2

MS of DAR1 ADC= The MS of G1/G2(Herceptin)+[(180.16)(Galactose)-18.02(H₂O)]+ [(1808.8)(GDP-Fuco-vc-PAB-MMAE)-443.20(GDP)]*1

MS of DAR0 ADC= The MS of G1/G2(Herceptin)+ [(180.16)(Galactose)-18.02(H₂O)]+[(1808.8)(GDP-Fuco-vc-PAB-MMAE)-443.20(GDP)]*0

2. Measure the DAR value of the conjugated antibody

DAR value=[Intensity (DAR0 ADC)*0+ Intensity (DAR1 ADC)*1+ Intensity (DAR2 ADC)*2+ Intensity (DAR3 ADC)*3+ Intensity (DAR4 ADC)*4] / [Intensity (DAR0 ADC)+ Intensity (DAR1 ADC)+ Intensity (DAR2 ADC)+ Intensity (DAR3 ADC)+ Intensity (DAR4 ADC)]

•Measure the DAR value of the conjugated antibody (HIC-HPLC based)

Column: TSKgel Butyl-NPR column (4.6 mm × 35 mm, 2.5 μm; TOSOH)

Mobile phase A: 20 mM sodium phosphate, 1.5 M ammonium sulfate (pH 6.9)

Mobile phase B: 75% (v/v) 20 mM sodium phosphate, 25% (v/v) isopropanol (pH 6.9)

Flow rate: 0.4 mL/min

Column temperature: 30 °C


Detection wavelength: 280nm

Loading Amount: 10 μg (5-10 μL)

Gradient:

Time	A (%)	B (%)
0.00	100	0
1.00	100	0
13.00	0	100
13.50	0	100
14.50	100	0
18.50	100	0

DAR value calculation:



DAR value = {[DAR0 (7.014) peak area percent]*0+[DAR1 (7.760) peak area percent]*1+[DAR2 (8.162) peak area percent]*2+[DAR3 (8.834) peak area percent]*3+[DAR4 (9.331) peak area percent]*4}/100

Trouble shooting:

Problem	Possible Cause	Solution
DAR value <3.2	The reaction does not reach the end point	Reduce the concentration of antibody, or increase the concentration of enzyme in the reaction system, or extend the reaction time.
	The antibody has some glycan types, which can not be modified, such as high mannose type	---
DAR value >4.0	The antibody has more than one N-linked glycosylation sites.	---
Precipitation occurred immediately when MnCl ₂ is added to the reaction system.	The buffer exchange of the antibody is incomplete.	Buffer exchange of the antibody again.
Precipitation occurred during the reaction.	A small amount of the enzyme inactivates during the reaction.	It has no effect on the final results.

Reference:

[1] Van Geel R, Wijdeven M A, Heesbeen R, Heesbeen R, *et al.* Chemoenzymatic conjugation of toxic payloads to the globally conserved N-glycan of native mAbs provides homogeneous and highly efficacious antibody-drug conjugates *J. Bioconjugate Chemistry*, 2015, 26(11):2233-2242.

[2] Yang Y, Zhentao S, Tian T, *et al.* Reducing the Complexity of Fc Glycan Enables the Construction of Antibody conjugates with Unexpected High Efficiency and Payload Capacity via Glycoengineering *Angewandte Chemie International Edition* doi: <https://doi.org/10.1101/2022.09.04.506510>.

Statement of Rights:

AGlink® products are co-developed by ACROBiosystems Co., Ltd. and Glyco-therapy Biotechnology Co., Ltd. based on Glyco-therapy's conjugation platform (YT Conju™) and are not authorized for commercial use. The sale and use of AGlink® products doesn't represent any forms of license or authorization of the attached technologies for developing any diagnostic, therapeutic or prophylactic products for commercial purposes.

Note of Application:

AGlink® products are intended for research use only by the Buyer or its users and can not be used for any other purposes, including, but not limited to, use for any clinical diagnostic type use, in vitro or in vivo therapeutic use, prophylactic use, any type of consumption or application to human or animals.