

## Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD32a-Medium Expression)

Catalog Number: CK-002

| Components   | Catalog No.   | Size   |
|--|---------------|--|
| Human 4-1BB (Luc) HEK293 Reporter Cell               | CHEK-ATF073   | 2 × (1 vial contains ~5×10 <sup>6</sup> cells) |
| CHO/Human CD32a Stable Cell Line (Medium Expression) | SCCHO-ATP061M | 2 × (1 vial contains ~5×10 <sup>6</sup> cells) |


### • Description

The Cell-based Kit consists of two engineered cell lines, Human 4-1BB (Luc) HEK293 Reporter Cell (Cat.No.CHEK-ATF073) and CHO/Human CD32a Stable Cell Line (Medium Expression) (Cat.No.SCCHO-ATP061M).

The Human 4-1BB (Luc) HEK293 Reporter Cell was engineered to not only express NF-κB signaling response element, but also express the receptor full length human 4-1BB (Gene ID: 3604), which can drive luciferase expressing systems by 4-1BB ligand/ agonist antibody stimulation. The CHO/Human CD32a Stable Cell Line was engineered to express full length human CD32a receptor (Gene ID: 2212), with different levels of CD32a expression (High, Medium, Low), which can be used to test agonist antibody whether in a CD32a-dependent manner to strengthen the agonistic activity. When co-cultured with Human 4-1BB (Luc) HEK293 Reporter Cell and anti-4-1BB agonist antibody, the anti-4-1BB antibody can be crosslinked, thereby strengthening 4-1BB pathway-activated luminescence.

### • Application

- Screen for Anti-human 4-1BB antibodies whether in a CD32a-dependent manner to strengthen the agonistic activity



## Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD32a-Medium Expression) Data Sheet

### • Cell Line Profile of Human 4-1BB (Luc) HEK293 Reporter Cell

|                               |   |
|-------------------------------|---|
| <b>Cell line</b>              | Human 4-1BB (Luc) HEK293 Reporter Cell      |
| <b>Host Cell</b>              | HEK293                                      |
| <b>Property</b>               | Adherent                                    |
| <b>Complete Growth Medium</b> | DMEM + 10% FBS                              |
| <b>Selection Marker</b>       | Puromycin (2 µg/mL) + Hygromycin (20 µg/mL) |
| <b>Incubation</b>             | 37°C with 5% CO <sub>2</sub>                |
| <b>Doubling Time</b>          | 22-24 hours                                 |
| <b>Transduction Technique</b> | Lentivirus                                  |

### • Cell Line Profile of CHO/Human CD32a Stable Cell Line (Medium Expression)

|                               |  |
|-------------------------------|--|
| <b>Cell line</b>              | CHO/Human CD32a Stable Cell Line (Medium Expression) |
| <b>Host Cell</b>              | CHO  |
| <b>Property</b>               | Adherent   |
| <b>Complete Growth Medium</b> | F-12K + 10% FBS                                      |
| <b>Selection Marker</b>       | Hygromycin (20 µg/mL)                                |
| <b>Incubation</b>             | 37°C with 5% CO <sub>2</sub>                         |
| <b>Doubling Time</b>          | 22-24 hours  |
| <b>Transduction Technique</b> | Lentivirus   |

• ***Cell Culture of Human 4-1BB (Luc) HEK293 Reporter Cell***

***Materials Required***

- DMEM medium (Gibco, Cat.No.11965-092)
- Fetal bovine serum (CellMax, Cat.No.SA211.02)
- Puromycin (InvivoGen, Cat.No.ant-pr-5b)
- Hygromycin B (Invitrogen, Cat.No.10687010)
- Complete Growth Medium: DMEM + 10% FBS
- Culture Medium: DMEM + 10% FBS, Hygromycin (20 µg/mL), Puromycin (2 µg/mL)
- Freeze Medium: 90% FBS, 10% (V/V) DMSO
- T-75 Culture flask (Corning, 430641)
- Cryogenic storage vials (SARSTEDT, 72.379.007)
- Thermostat water bath
- Centrifuge
- Luna cell counter (Logos Biosystems, LUNA- II )
- CO<sub>2</sub> Incubator (Thermo, 3111)
- Biological Safety Cabinet (Thermo, 1389)

***Recovery***

1. Thaw the vial by gentle agitation in a 37°C water bath. To reduce the possibility of contamination, keep the cap out of the water. Thawing should be rapid (approximately 2 minutes).
2. Remove the vial from the water bath as soon as the contents are thawed, and decontaminate by spraying with 70% ethanol. All the operations from this point on should be carried out under strict aseptic conditions.
3. Transfer the vial contents to a centrifuge tube containing 4.0 mL complete growth medium and spin at approximately 1000 rpm for 5 minutes.
4. Resuspend cell pellet with 5 mL complete growth medium and transfer the cell suspension into T-75 flask containing 10-15 mL of pre-warmed complete growth medium.
5. Incubate at 37°C with 5% CO<sub>2</sub> incubator until the cells are ready to be split.

### *Subculture*

1. Remove and discard culture medium.
2. Wash the cells once with sterile PBS.
3. Add 3 mL of 0.25% trypsin to cell culture flask. Place the flask at 37°C for 5-7 minutes, until 90% of the cells have detached.
4. Add 6.0 to 8.0 mL of culture medium and aspirate cells by gently pipetting.
5. Add appropriate aliquots of the cell suspension to new culture vessel.
6. Incubate at 37°C with 5% CO<sub>2</sub> incubator.

**Subcultivation Ratio:** A subcultivation ratio of 1:6 to 1:10 is recommended.

**Medium Renewal:** Every 2 to 3 days.

### *Cryopreservation*

1. Remove and discard spent medium.
2. Detach cells from the cell culture flasks with 0.25% trypsin.
3. Centrifuge at 1000 rpm for 5 min at RT to pellet cells.
4. Resuspend the cell pellets with culture medium and count viable cells.
5. Centrifuge at 1000 rpm for 5 min at RT and resuspend cells in freezing medium to a concentration of  $5 \times 10^6$  to  $1 \times 10^7$  cells/mL.
6. Aliquot into cryogenic storage vials. Place vials in a programmable cooler or an insulated box placed in a -80°C freezer overnight, then transferring to liquid nitrogen storage.

### *Storage*

- **Product format:** Frozen
- **Storage conditions:** Liquid nitrogen immediately upon receipt

• ***Cell Culture of CHO/Human CD32a Stable Cell Line (Medium Expression)***

***Materials Required***

- F-12K Nutrient Mixture (Gibco, Cat.No.21127-022)
- Fetal bovine serum (CellMax, Cat.No.SA211.02)
- Hygromycin B (Invitrogen, Cat.No.10687010)
- Complete Growth Medium: F-12K + 10% FBS
- Culture Medium: F-12K + 10% FBS, Hygromycin (20 µg/mL)
- Freeze Medium: 90% FBS, 10% (V/V) DMSO
- T-75 Culture flask (Corning, 430641)
- Cryogenic storage vials (SARSTEDT, 72.379.007)
- Thermostat water bath
- Centrifuge
- Luna cell counter (Logos Biosystems, LUNA- II )
- CO<sub>2</sub> Incubator (Thermo, 3111)
- Biological Safety Cabinet (Thermo, 1389)

***Recovery***

1. Thaw the vial by gentle agitation in a 37°C water bath. To reduce the possibility of contamination, keep the cap out of the water. Thawing should be rapid (approximately 2 minutes).
2. Remove the vial from the water bath as soon as the contents are thawed, and decontaminate by spraying with 70% ethanol. All the operations from this point on should be carried out under strict aseptic conditions.
3. Transfer the vial contents to a centrifuge tube containing 4.0 mL complete growth medium and spin at approximately 1000 rpm for 5 minutes.
4. Resuspend cell pellet with 5 mL complete growth medium and transfer the cell suspension into T-75 flask containing 10-15 mL of pre-warmed complete growth medium.
5. Incubate at 37°C with 5% CO<sub>2</sub> incubator until the cells are ready to be split.

### *Subculture*

1. Remove and discard culture medium.
2. Wash the cells once with sterile PBS.
3. Add 3 mL of 0.25% trypsin to cell culture flask. Place the flask at 37°C for 5-7 minutes, until 90% of the cells have detached.
4. Add 6.0 to 8.0 mL of culture medium and aspirate cells by gently pipetting.
5. Add appropriate aliquots of the cell suspension to new culture vessel.
6. Incubate at 37°C with 5% CO<sub>2</sub> incubator.

**Subcultivation Ratio:** A subcultivation ratio of 1:6 to 1:10 is recommended.

**Medium Renewal:** Every 2 to 3 days.


### *Cryopreservation*

1. Remove and discard spent medium.
2. Detach cells from the cell culture flasks with 0.25% trypsin.
3. Centrifuge at 1000 rpm for 5 min at RT to pellet cells.
4. Resuspend the cell pellets with culture medium and count viable cells.
5. Centrifuge at 1000 rpm for 5 min at RT and resuspend cells in freezing medium to a concentration of  $5 \times 10^6$  to  $1 \times 10^7$  cells/mL.
6. Aliquot into cryogenic storage vials. Place vials in a programmable cooler or an insulated box placed in a -80°C freezer overnight, then transferring to liquid nitrogen storage.


### *Storage*

- **Product format:** Frozen
- **Storage conditions:** Liquid nitrogen immediately upon receipt

• *Receptor Assay*



**Fig1. Expression analysis of human 4-1BB on Human 4-1BB (Luc) HEK293 Reporter Cell by FACS.** Cell surface staining was performed on Human 4-1BB (Luc) HEK293 Reporter Cell or negative control cell using PE-labeled anti-human 4-1BB antibody.




| Catalog No.   | Stable Cell Line                                     | MFI for CD32a (PE) |
|---------------|--|--------------------|
| SCCHO-ATP061L | CHO/Human CD32a Stable Cell Line (Low Expression)    | 5196.88            |
| SCCHO-ATP061M | CHO/Human CD32a Stable Cell Line (Medium Expression) | 11336.11           |
| SCCHO-ATP061H | CHO/Human CD32a Stable Cell Line (High Expression)   | 18309.90           |

**Fig2. Expression analysis of human CD32a on CHO/Human CD32a Stable Cell Line by FACS.** Cell surface staining using PE-labeled anti-human CD32a antibody was performed on CHO/Human CD32a Stable Cell Line with different expression levels: CHO/Human CD32a Stable Cell Line (Low Expression); CHO/Human CD32a Stable Cell Line (Medium Expression); CHO/Human CD32a Stable Cell Line (High Expression).



• *Application*

**CHO/CD32a (Medium Expression) Crosslinking**



**Fig3. Bioactivity analysis of anti-human 4-1BB antibody through CHO/Human CD32a Stable Cell Line (Medium Expression) crosslinking to test whether in a CD32a-dependent manner to strengthen the agonistic activity. The EC50 of anti-human 4-1BB antibody is approximately 0.05 µg/mL through CHO/Human CD32a Stable Cell Line (Medium Expression) crosslinking.**

## Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD32a-Medium Expression) Data Sheet

### • *License Disclosure*

This reporter cell is provided for research use only. This license does not permit you to share, distribute, sell, sublicense, or otherwise make this reporter cell available for use to other laboratories, departments, research institutions, hospitals, universities, or biotech companies. The license does not permit modification of this reporter cell in any way. Inappropriate use or distribution of this reporter cell will result in revocation of the license. Modifications of this cell line, transfer to another facility, or commercial use of the cells may require a separate license and additional fees. AcroBiosystems does not warrant the suitability of this reporter cell for any particular use, and does not accept any liability in connection with the handling or use of this reporter cell.

### • *Related Products*

| <u>Products</u>   | <u>Cat.No.</u> |
|---|----------------|
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD32a-Low Expression)           | CK-001         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD32a-High Expression)          | CK-003         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD16a (158V)-Low Expression)    | CK-004         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD16a (158V)-Medium Expression) | CK-005         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD16a (158V)-High Expression)   | CK-006         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD32b-Low Expression)           | CK-007         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD32b-Medium Expression)        | CK-008         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD32b-High Expression)          | CK-009         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD64-Low Expression)            | CK-010         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD64-Medium Expression)         | CK-011         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human CD64-High Expression)           | CK-012         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human PD-L1-Low Expression)           | CK-013         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human PD-L1-Medium Expression)        | CK-014         |
| Cell-based Screening Kit for Anti-human 4-1BB Antibody (Human PD-L1-High Expression)          | CK-015         |