



Synonym

PVR,FLJ25946,PVS,CD155,TAGE4,HVED,NECL5

Source

Human CD155, Fc Tag(CD5-H5251) is expressed from human 293 cells (HEK293). It contains AA Trp 21 - Asn 343 (Accession # [NP_006496.2](#)).

Predicted N-terminus: Trp 21

Molecular Characterization

CD155(Trp 21 - Asn 343) NP_006496.2	Fc(Pro 100 - Lys 330) P01857
--	---------------------------------

This protein carries a human IgG1 Fc tag at the C-terminus.

The protein has a calculated MW of 61.7 kDa. The protein migrates as 70-90 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>95% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in

Tris with Glycine, Arginine and NaCl, pH7.5 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE




Human CD155, Fc Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95% (With [Star Ribbon Pre-stained Protein Marker](#)).


Bioactivity-ELISA

Discounts, Gifts,
and more!






Immobilized Human TIGIT, His Tag (Cat. No. TIT-H52H3) at 2 µg/mL (100 µL/well) can bind Human CD155, Fc Tag (Cat. No. CD5-H5251) with a linear range of 8-128 ng/mL (QC tested).




Serial dilutions of Human TIGIT Neutralizing antibody were added into Human CD155, Fc Tag (Cat. No. CD5-H5251): Biotinylated Human TIGIT, Fc,Avitag (Cat. No. TIT-H82F1) binding reactions. The half maximal inhibitory concentration (IC50) is 0.06065 µg/mL (Routinely tested).

Bioactivity-BLI





Loaded Human CD155, Fc Tag (Cat. No. CD5-H5251) on Protein A Biosensor, can bind Human TIGIT, His Tag with an affinity constant of 0.23 µM as determined in BLI assay (ForteBio Octet Red96e) (Routinely tested).



Loaded Human CD155, Fc Tag (Cat. No. CD5-H5251) on Protein A Biosensor, can bind Cynomolgus DNAM-1, His Tag (Cat. No. DN1-C52H9) with an affinity constant of 0.88 µM as determined in BLI assay (ForteBio Octet Red96e) (Routinely tested).

Bioactivity-FACS



Discounts, Gifts, and more!



Human CD155 / PVR Protein, Fc Tag

Catalog # CD5-H5251



BIOSYSTEMS
Acro

FACS assay shows that Human CD155, Fc Tag (Cat. No. CD5-H5251) can bind to 293T cell overexpressing human TIGIT. The concentration of CD155 is 1 µg/mL (Routinely tested).

FACS analysis shows that the binding of Human CD155, Fc Tag (Cat. No. CD5-H5251) to 293T overexpressing TIGIT was inhibited by increasing concentration of neutralizing Anti-Human TIGIT MAbs. The concentration of CD155 used is 1 µg/mL. The IC₅₀ is 0.09422 µg/mL (Routinely tested).

Background

CD155 is a Type I transmembrane glycoprotein in the immunoglobulin superfamily. Commonly known as Poliovirus Receptor (PVR) due to its involvement in the cellular poliovirus infection in primates, CD155's normal cellular function is in the establishment of intercellular adherens junctions between epithelial cells.

CD155/PVR was originally isolated based on its ability to mediate polio virus attachment to host cells. The fulllength (or CD155 alpha isoform) is synthesized as a 417 amino acid (aa) precursor that contains a 20 aa signal sequence, a 323 aa extracellular region, a 24 aa TM segment and a 50 aa cytoplasmic tail. The extracellular region contains one N terminal V type and two C2 type Ig like domains.

CD155 is a transmembrane protein with 3 extracellular immunoglobulin-like domains, D1-D3, where D1 is recognized by the virus. Low resolution structures of CD155 complexed with poliovirus have been obtained using electron microscopy while a high resolution structures of the ectodomain D1 and D2 of CD155 were solved by x-ray crystallography.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

Discounts, Gifts,
and more!



➤ www.acrobiosystems.com

12/20/2023