Catalog # APL-C52D2



Synonym

APRIL, TNFSF13, TALL-2, TRDL-1, CD256, TALL2, ZTNF2

Source

Cynomolgus APRIL, His,Flag Tag (APL-C52D2) is expressed from human 293 cells (HEK293). It contains AA Gln 111 - Leu 250 (Accession # <u>XP_005582818.1</u>) trimer Design. Predicted N-terminus: His

Molecular Characterization

This protein carries a polyhistidine tag at the N-terminus, followed by a flag tag

The protein has a calculated MW of 50.4 kDa. The protein migrates as 55-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per μg by the LAL method.

Purity

>95% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Cynomolgus APRIL, His,Flag Tag on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-MALS



The purity of Cynomolgus APRIL, His,Flag Tag (Cat. No. APL-C52D2) is more than 90% and the molecular weight of this protein is around 57-69 kDa verified by SEC-MALS.

<u>Report</u>

Bioactivity-ELISA


>>> www.acrobiosystems.com


3/15/2023

Cynomolgus APRIL / TNFSF13 Protein, His, Flag Tag, active trimer (MALS verified)



Catalog # APL-C52D2





Immobilized Cynomolgus APRIL, His,Flag Tag (Cat. No. APL-C52D2) at 1 μ g/mL (100 μ L/well) can bind Human TACI (R119A, R122A), Fc Tag (Cat. No. TAI-H5253) with a linear range of 2-13 ng/mL (QC tested).

Immobilized Cynomolgus APRIL, His,Flag Tag (Cat. No. APL-C52D2) at 1 μ g/mL (100 μ L/well) can bind Human BCMA, Fc Tag (HPLC verified) (Cat. No. BC7-H5254) with a linear range of 0.1-1 ng/mL (Routinely tested).

Background

APRIL(a proliferation-inducing ligand) is also known as Tumor necrosis factor ligand superfamily member 13, TALL-2, TRDL-1, CD256, TNFFSF 13, cytokine that binds to TNFRSF13B/TACI and to TNFRSF17/BCMA. APRIL is a cytokine of the tumor necrosis factor family associated mainly with hematologic malignancies. The closely related TNF family ligands B cell activation factor (BAFF) and a proliferation-inducing ligand (APRIL) serve in the generation and maintenance of mature B-lymphocytes. Both BAFF and APRIL assemble as homotrimers that bind and activate several receptors that they partially share. BAFF-APRIL heteromers of different stoichiometries have distinct receptor-binding properties and activities. In addition, expression of APRIL was regulated by miR-145 in GC cells.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.



3/15/2023