Catalog # CD4-H82E9



Synonym

CD24,CD24A

Source

Biotinylated Human CD24, His,Avitag(CD4-H82E9) is expressed from human 293 cells (HEK293). It contains AA Ser 27 - Gly 59 (Accession # <u>P25063-1</u>). Predicted N-terminus: Ser 27

Molecular Characterization

CD24(Ser 27 - Gly 59) P25063-1 Poly-his Avi

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (AvitagTM).

The protein has a calculated MW of 6.7 kDa. The protein migrates as 18-33 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

Biotinylation of this product is performed using Avitag[™] technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Endotoxin

Less than 1.0 EU per μg by the LAL method.

SDS-PAGE

kDa	М	R
116.0		
66.2	—	
45.0	-	
35.0	-	-
25.0	_	
18.4	-	
14.4	-	

Biotinylated Human CD24, His, Avitag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

Bioactivity-ELISA


>>> www.acrobiosystems.com

10/8/2023

Biotinylated Human CD24 Protein, His,Avitag™



Catalog # CD4-H82E9



Immobilized Biotinylated Human CD24, His, Avitag (Cat. No. CD4-H82E9) at 5 μ g/mL (100 μ L/well) on streptavidin precoated (0.5 μ g/well) plate, can bind Anti-CD24 MAb (SN3) with a linear range of 0.6-10 ng/mL (QC tested).

Background

CD24 may have a pivotal role in cell differentiation of different cell types. Signaling could be triggered by the binding of a lectin-like ligand to the CD24 carbohydrates, and transduced by the release of second messengers derived from the GPI-anchor. Modulates B-cell activation responses. Promotes AG-dependent proliferation of B-cells, and prevents their terminal differentiation into antibody-forming cells. In association with SIGLEC10 may be involved in the selective suppression of the immune response to danger-associated molecular patterns (DAMPs) such as HMGB1, HSP70 and HSP90. Plays a role in the control of autoimmunity.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

