

**Synonym**

CLEC12A,MICL,CLL-1,CLL1,DCAL2,DCAL-2,CD371

Source

Biotinylated Human CLEC12A, Fc,Avitag(CLA-H82F7) is expressed from human 293 cells (HEK293). It contains AA His 65 - Ala 265 (Accession # [Q5QGZ9-2](#)).

Predicted N-terminus: Gly

Molecular Characterization

Avi	Fc(Pro 100 - Lys 330) P01857	CLEC12A(His 65 - Ala 265) Q5QGZ9-2
-----	---------------------------------	---------------------------------------

This protein carries an Avi tag (Avitag™) at the N-terminus, followed by a human IgG1 Fc tag.

The protein has a calculated MW of 51.9 kDa. The protein migrates as 65-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

Biotinylation of this product is performed using Avitag™ technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>95% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.


Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.



This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE

Biotinylated Human CLEC12A, Fc,Avitag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.


Bioactivity-ELISA**Discounts, Gifts,
and more!**



Immobilized Biotinylated Human CLEC12A, Fc,Avitag (Cat. No. CLA-H82F7) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Monoclonal Anti-Human CLEC12A Antibody, Human IgG1 with a linear range of 0.1-2 ng/mL (QC tested).

Immobilized Monoclonal Anti-Human CLEC12A Antibody, Human IgG1 at 1 µg/mL (100 µL/well) can bind Biotinylated Human CLEC12A, Fc,Avitag (Cat. No. CLA-H82F7) with a linear range of 0.8-13 ng/mL (Routinely tested).

Bioactivity-FACS



2e5 of anti-CLEC12A CAR-293 were stained with 100 µL of 3 µg/mL of Biotinylated Human CLEC12A, Fc,Avitag (Cat. No. CLA-H82F7) and negative control protein respectively, washed and then followed by PE-SA and analyzed with FACS (Routinely tested).

Background

CLEC12A (C-type lectin domain family 12 member A) is also known as CLL1, DCAL2, MICL. Clec12a is an inhibitory receptor for uric acid crystals that regulates inflammation in response to cell death. Cell surface receptor that modulates signaling cascades and mediates tyrosine phosphorylation of target MAP kinases. Evidence of distinct disease propagating stem cells in myelodysplastic syndrome (MDS) has emerged in recent years. The role of CLEC12A in MDS, however, remains to be elucidated. Furthermore, CLEC12A has been proposed as a promising marker of leukaemic stem cells in AML.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

Discounts, Gifts, and more!

