

Synonym

Gn protein, glycoprotein polyprotein, Phlebovirus_G1, Phlebovirus_G2

Source

SFTS virus Gn Protein (Human/China/HB29/2010), His Tag(GNN-S52H3) is expressed from human 293 cells (HEK293). It contains AA Asp 20 - Lys 452 (Accession # [F1BA47-1](#)).

Predicted N-terminus: Asp 20

Molecular Characterization

Gn(Asp 20 - Lys 452)
F1BA47-1 Poly-his

This protein carries a polyhistidine tag at the C-terminus

The protein has a calculated MW of 49.1 kDa. The protein migrates as 50-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>95% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.


Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.


Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE

SFTS virus Gn Protein (Human/China/HB29/2010), His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-MALS

The purity of SFTS virus Gn Protein (Human/China/HB29/2010), His Tag (Cat. No. GNN-S52H3) is more than 85% and the molecular weight of this protein is around 95-115 kDa verified by SEC-MALS.

[Report](#)

Background

Severe fever with thrombocytopenia syndrome (SFTS) is an emerging viral hemorrhagic fever (VHF) endemic to China, South Korea, Japan, and Vietnam. Severe fever with thrombocytopenia syndrome (SFTS) is an infectious disease with a high fatality rate, caused by SFTS virus (SFTSV). To our knowledge, no efficient SFTSV vaccine exists.

Clinical and Translational Updates

SFTS virus Gn Protein (Human/China/HB29/2010), His Tag (MALS verified)

Catalog # GNN-S52H3



Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.