



Synonym

HLA-A*0201 & B2M & HPV16-E6

Source

PE-Labeled Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Tetramer Protein(HLH-HP2H7) is expressed from human 293 cells (HEK293). It contains AA Gly 25 - Ile 308 (HLA-A*02:01) & Ile 21 - Met 119 (B2M) & TIHDIILECV peptide (Accession # [AAA59606.1](#) (HLA-A*02:01) & [P61769-1](#) (B2M) & TIHDIILECV).

Predicted N-terminus: Gly 25 & Ile 21

Molecular Characterization

PE-Labeled Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Tetramer Protein is assembled by biotinylated monomer and PE-labeled streptavidin.

Biotinylated Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Complex Protein is produced by co-expression of HLA and B2M loaded with HPV16-E6 peptide. Biotinylated Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Complex Protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag™).

Conjugate

PE

Excitation Wavelength: 488 nm / 561 nm

Emission Wavelength: 575 nm

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Supplied as 0.2 µm filtered solution in PBS, 1% BSA, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Shipping

This product is supplied and shipped with dry ice, please inquire the shipping cost.

Storage


Please protect from light and avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- The product MUST be stored at -70°C or lower upon receipt;
- -70°C for 3 months under sterile conditions.

Bioactivity-ELISA

PE-Labeled Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Tetramer Protein ELISA
0.1 µg of PE-Labeled Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Tetramer Protein per well



Immobilized PE-Labeled Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Tetramer Protein (Cat. No. HLH-HP2H7) at 1 µg/mL (100

Discounts, Gifts, and more!



PE-Labeled Human HLA-A*02:01&B2M&HPV16-E6 (TIHDIILECV) Tetramer Protein

Catalog # HLH-HP2H7



BIOSYSTEMS
Acro

μL/well) can bind Anti-B2M Antibody, Human IgG1 with a linear range of 0.1-2 ng/mL (Routinely tested).

Background

Human papillomavirus (HPV) is A kind of papillomavirus belonging to the milk polypoid virus family. It is a spherical D virus, which can cause the proliferation of squamous epithelium of human skin mucosa. HPV(human papillomavirus) for common warts, genital warts (condyloma acuminatum), and other symptoms. There are many types of human papillomavirus (HPV), with HPV 16 and 18 being high-risk types known to significantly increase the risk of cervical, vaginal and vulvar cancers in women and men. The PE-labeled human HLA-A*0201 HPV (TIHDIILECV) tetramer protein is a complex of HLA-A*0201 of the MHC Class I, B2M, and TIHDIILECV peptide of the HPV.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

Discounts, Gifts,
and more!



➤ www.acrobiosystems.com

11/15/2023