



## Synonym

KREMEN2, Kremen-2, KRM2, Dickkopf receptor 2, Kringle domain-containing transmembrane protein 2

## Source

Human Kremen-2, His Tag(KR2-H52H3) is expressed from human 293 cells (HEK293). It contains AA Gly 26 - Ala 364 (Accession # [Q8NCW0-1](#)).

## Molecular Characterization

KRM2(Gly 26 - Ala 364)  
Q8NCW0-1 Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 38.0 kDa. The protein migrates as 45-60 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

## Endotoxin

Less than 1.0 EU per µg by the LAL method.

## Purity

>95% as determined by SDS-PAGE.

## Formulation

Lyophilized from 0.22 µm filtered solution in PBS, 0.2M Arginine, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

## Reconstitution

Please see Certificate of Analysis for specific instructions.

*For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.*

## Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

*Please avoid repeated freeze-thaw cycles.*

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

## SDS-PAGE




Human Kremen-2, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95% (With [Star Ribbon Pre-stained Protein Marker](#)).

## Bioactivity-ELISA

Discounts, Gifts,  
and more!





Immobilized Human Dkk-1, Fc Tag (Cat. No. DK1-H5258) at 2 µg/mL (100 µL/well) can bind Human Kremen-2, His Tag (Cat. No. KR2-H52H3) with a linear range of 10-78 ng/mL (QC tested).

## Background

Kremen2 (Krm2) plays an important role in embryonic development, bone formation, and tumorigenesis as a crucial regulator of classical Wnt/ $\beta$ -catenin signaling pathway. Krm1 and its paralog Krm2 share the ability to bind Dkk1 and inhibit Wnt signaling, both processes previously shown to rely on the extracellular domain. Previous studies have shown that Krm2 may be involved in embryonic development, bone formation, neural ridge formation and tumorigenesis and could be a biomarker of grading and a potential therapeutic target in gastric cancer.

## Clinical and Translational Updates

Discounts, Gifts,  
and more!

