

#### Synonym

LILRB4,ILT3,LIR5,CD85K,HM18

#### Source

Human LILRB4, Fc Tag(LI4-H5256) is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Glu 259 (Accession # <u>AAH26309.1</u>). Predicted N-terminus: Gln 22

#### **Molecular Characterization**

LILRB4(Gln 22 - Glu 259) Fc(Pro 100 - Lys 330)
AAH26309.1 P01857

This protein carries a human IgG1 Fc tag at the C-terminus.

The protein has a calculated MW of 52.7 kDa. The protein migrates as 60-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

### Endotoxin

Less than 1.0 EU per µg by the LAL method.

## **Purity**

>95% as determined by SDS-PAGE.

>95% as determined by SEC-MALS.

#### **Formulation**

Lyophilized from  $0.22~\mu m$  filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

#### Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

### **Storage**


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


# SDS-PAGE



Human LILRB4, Fc Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.

## **Bioactivity-ELISA**


### **SEC-MALS**




The purity of Human LILRB4, Fc Tag (Cat. No. LI4-H5256) is more than 95% and the molecular weight of this protein is around 103-123 kDa verified by SEC-MALS.

Report

## Human LILRB4 / CD85k / ILT3 Protein, Fc Tag (MALS verified)







Immobilized Human LILRB4, Fc Tag (Cat. No. LI4-H5256) at 1  $\mu$ g/mL (100  $\mu$ L/well) can bind Anti-LILRB4 Antibody, Human IgG1 with a linear range of 0.4-6 ng/mL (QC tested).

## **Background**

Leukocyte immunoglobulin-like receptor subfamily B member 4 (LILRB4) is also known as CD85 antigen-like family member K (CD85K), Immunoglobulin-like transcript 3 (ILT-3), Leukocyte immunoglobulin-like receptor 5 (LIR-5), Monocyte inhibitory receptor HM18, which belongs to the leukocyte immunoglobulin-like receptor (LIR) family. LILRB4 / CD85K contains 2 Ig-like C2-type (immunoglobulin-like) domains. CD85K is detected in monocytes, macrophages, dendritic cells, lung, natural killer cells and B-cells. LILRB4 / CD85K is receptor for class I MHC antigens. CD85K recognizes a broad spectrum of HLA-A, HLA-B, HLA-C and HLA-G alleles, involved in the down-regulation of the immune response and the development of tolerance. LILRB4 interferes with TNFRSF5-signaling and NF-kappa-B up-regulation and inhibits receptor-mediated phosphorylation of cellular proteins and mobilization of intracellular calcium ions.

#### **Clinical and Translational Updates**

Please contact us via <u>TechSupport@acrobiosystems.com</u> if you have any question on this product.