

Synonym

NOTCH2,hN2,N2ECD

Source

Human NOTCH2, Fc Tag(NO2-H5255) is expressed from human 293 cells (HEK293). It contains AA Leu 26 - Gln 530 (Accession # [Q04721-1](#)).

Predicted N-terminus: Leu 26

Molecular Characterization

NOTCH2(Leu 26 - Gln 530) Q04721-1	Fc(Pro 100 - Lys 330) P01857
--------------------------------------	---------------------------------

This protein carries a human IgG1 Fc tag at the C-terminus

The protein has a calculated MW of 80.5 kDa. The protein migrates as 90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method.

Purity

>95% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.


Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.


Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


SDS-PAGE

Human NOTCH2, Fc Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.


Bioactivity-ELISA**SEC-MALS**

The purity of Human NOTCH2, Fc Tag (Cat. No. NO2-H5255) is more than 85% and the molecular weight of this protein is around 175-190 kDa verified by SEC-MALS.

[Report](#)




Immobilized Human DLL4, His Tag (Cat. No. DL4-H5227) at 5 µg/mL (100 µL/well) can bind Human NOTCH2, Fc Tag (Cat. No. NO2-H5255) with a linear range of 5-78 ng/mL (QC tested).



Immobilized Human NOTCH2, Fc Tag (Cat. No. NO2-H5255) at 1 µg/mL (100 µL/well) can bind Bispecific Antibody (NOTCH2 x NOTCH3) with a linear range of 1-20 ng/mL (Routinely tested).

Bioactivity-BLI



Loaded Human NOTCH2, Fc Tag (Cat. No. NO2-H5255) on Protein A Biosensor, can bind Human DLL4, His Tag (Cat. No. DL4-H5227) with an affinity constant of 77.6 nM as determined in BLI assay (ForteBio Octet Red96e) (Routinely tested).

Background

Functions as a receptor for membrane-bound ligands Jagged-1 (JAG1), Jagged-2 (JAG2) and Delta-1 (DLL1) to regulate cell-fate determination. Upon ligand activation through the released notch intracellular domain (NICD) it forms a transcriptional activator complex with RBPJ/RBPSUH and activates genes of the enhancer of split locus. Affects the implementation of differentiation, proliferation and apoptotic programs (By similarity). Involved in bone remodeling and homeostasis. In collaboration with RELA/p65 enhances NFATc1 promoter activity and positively regulates RANKL-induced osteoclast differentiation. Positively regulates self-renewal of liver cancer cells.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.