

Synonym

Catalog # PDL-C82E8

PD-L1,CD274,B7-H1,PDCD1L1,PDCD1LG1

Source

Biotinylated Cynomolgus / Rhesus macaque PD-L1, His,Avitag(PDL-C82E8) is expressed from human 293 cells (HEK293). It contains AA Phe 19 - Arg 238 (Accession # <u>G7PSE7-1</u>).

Predicted N-terminus: Phe 19

Molecular Characterization



This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (AvitagTM).

The protein has a calculated MW of 28.8 kDa. The protein migrates as 32-40 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

Biotinylation of this product is performed using Avitag[™] technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.


Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Endotoxin

Less than 1.0 EU per μ g by the LAL method.

SDS-PAGE



Biotinylated Cynomolgus / Rhesus macaque PD-L1, His, Avitag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The

Purity

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 μ m filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.





The purity of Biotinylated Cynomolgus / Rhesus macaque PD-L1, His,Avitag (Cat. No. PDL-C82E8) is more than 90% and the molecular weight of this

purity of the protein is greater than 90%.

Bioactivity-ELISA


protein is around 31-47 kDa verified by SEC-MALS. Report

>> www.acrobiosystems.com


9/22/2023




Catalog # PDL-C82E8




Immobilized Biotinylated Cynomolgus / Rhesus macaque PD-L1, His, Avitag (Cat. No. PDL-C82E8) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Human PD-1, Fc Tag, low endotoxin (Cat. No. PD1-H5257) with a linear range of 0.3-2 ng/mL (QC tested).



Immobilized Biotinylated Cynomolgus / Rhesus macaque PD-L1, His, Avitag (Cat. No. PDL-C82E8) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Cynomolgus PD-1, Fc Tag (Cat. No. PD1-C5254) with a linear range of 0.3-5 ng/mL (Routinely tested).



Immobilized Biotinylated Cynomolgus / Rhesus macaque PD-L1, His, Avitag (Cat. No. PDL-C82E8) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Anti-PD-L1 Neutralizing Antibody with a linear range of 0.1-3 ng/mL (Routinely tested).





Biotinylated Cynomolgus / Rhesus macaque PD-L1, His, Avitag ELISA

Biotinylated Cynomolgus / Rhesus macaque PD-LI, His,Avitag Conc. (ng/mL)

Immobilized Anti-PD-L1 Neutralizing Antibody at 1 µg/mL (100 µL/well) can bind Biotinylated Cynomolgus / Rhesus macaque PD-L1, His, Avitag (Cat. No. PDL-C82E8) with a linear range of 0.2-10 ng/mL (Routinely tested).


Bioactivity-SPR

>> www.acrobiosystems.com





Catalog # PDL-C82E8



Cynomolgus PD-1, Fc Tag (Cat. No. PD1-C5254) captured on CM5 chip via Anti-human IgG Fc antibodies surface can bind Biotinylated Cynomolgus / Rhesus macaque PD-L1, His,Avitag (Cat. No. PDL-C82E8) with an affinity constant of 6.59 μ M as determined in a SPR assay (Biacore 8K) (Routinely tested).

Background

Programmed cell death 1 ligand 1 (PDL1) is also known as B7-H, B7H1, MGC142294, MGC142296, PD-L1, PDCD1L1 and PDCD1LG1, which is a member of the growing B7 family of immune molecules and is involved in the regulation of cellular and humoral immune responses. PDL1 is a cell surface immunoglobulin superfamily with two Ig-like domains within the extracellular region and a short cytoplasmic domain. This protein is broadly expressed in the majority of peripheral tissues as well as hematopoietic cells. Interaction between PDL1 and its receptors belonging to the CD28 family of molecules provide both stimulatory and inhibitory signals in regulating T cell activation and tolerance. PDL1 may inhibit ongoing T-cell responses by inducing apoptosis and arresting cell-cycle progression.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.



