

Source

Monoclonal Anti-RSV-Pre-F0 specific Antibody, Human IgG1 (12C6) is a chimeric monoclonal antibody recombinantly expressed from CHO cells, which combines the variable region of a mouse monoclonal antibody with human IgG1 constant domain. The mouse monoclonal antibody is produced from a hybridoma resulting from fusion of SP2/0 myeloma and B-lymphocytes obtained from a mouse immunized with RSV-Pre-F0.

Isotype

Human IgG1 | Human Kappa

Specificity

This product is a specific antibody specifically reacts with RSV-Pre-F0.

Application

ELISA

Formulation

Please contact us for detailed information.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.


Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70 °C for 3 months after reconstitution;


Bioactivity-Elisa

Monoclonal Anti-RSV-Pre-F0 specific Antibody, Human IgG1 (12C6) ELISA
0.1 µg of HRSV (A) Pre-fusion glycoprotein F0, His Tag per well



Immobilized HRSV (A) Pre-fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H7) at 1 µg/mL (100 µL/well) can bind Monoclonal Anti-RSV-Pre-F0 specific Antibody, Human IgG1 (12C6) (Cat. No. RS0-S286) with a linear range of 0.2-2 ng/mL. HRSV (A) Post-fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H6) is verified not recognized by Monoclonal Anti-RSV-Pre-F0 specific Antibody, Human IgG1 (12C6) (Cat. No. RS0-S286) in low concentration (QC tested).

Monoclonal Anti-RSV-Pre-F0 specific Antibody, Human IgG1 (12C6) ELISA
0.1 µg of HRSV (A) Pre-fusion glycoprotein F0, His Tag per well



Immobilized HRSV (A) Pre-fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H7) at 1 µg/mL (100 µL/well) can bind Monoclonal Anti-RSV-Pre-F0 specific Antibody, Human IgG1 (12C6) (Cat. No. RS0-S286)/Anti-Fusion glycoprotein F0 Antibody, Human IgG1 (D25) with a linear range of 0.1-2 ng/mL. HRSV (A) Post-fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H6) is verified not recognized by Monoclonal Anti-RSV-Pre-F0 specific Antibody, Human IgG1 (12C6) (Cat. No. RS0-S286)/Anti-Fusion glycoprotein F0 Antibody, Human IgG1 (D25) in low concentration (Routinely tested).

Background

Monoclonal Anti-RSV-Pre-F0 specific Antibody, Human IgG1 (12C6)

Catalog # RS0-S286



Human respiratory syncytial virus (HRSV) is the most common etiological agent of acute lower respiratory tract disease in infants and can cause repeated infections throughout life. The RSV fusion glycoprotein (RSV F) is the principal target of RSV neutralizing antibodies in human sera. The RSV F is a type I viral fusion protein synthesized as inactive, single-chain polypeptides that assemble into trimers. RSV F fuses the viral and host cell membranes by irreversible protein refolding from the labile prefusion conformation to the stable post-fusion conformation.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.