

HRSV (A) Pre-fusion glycoprotein F0, His Tag (MALS verified)

Catalog # RSF-V52H7



Source

HRSV (A) Fusion glycoprotein F0, His Tag (RSF-V52H7) is expressed from human 293 cells (HEK293).

Predicted N-terminus: Gln 26

Molecular Characterization

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 56.1 kDa. The protein migrates as 60-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per μg by the LAL method.

Purity

>95% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Supplied as 0.2 μm filtered solution in PBS with trehalose as protectant.

Contact us for customized product form or formulation.

Shipping

This product is supplied and shipped with dry ice, please inquire the shipping cost.


Storage

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- The product MUST be stored at -70°C or lower upon receipt;
- -70°C for 3 months under sterile conditions.


SDS-PAGE



HRSV (A) Pre-fusion glycoprotein F0, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.

Bioactivity-ELISA

SEC-MALS



The purity of HRSV (A) Pre-fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H7) is more than 90% and the molecular weight of this protein is around 175-200 kDa verified by SEC-MALS.

[Report](#)

Discounts, Gifts,
and more!




HRSV (A) Pre-fusion glycoprotein F0, His Tag (MALS verified)


Catalog # RSF-V52H7




BIOSYSTEMS
Acro




Immobilized HRSV (A) Fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H7) at 2 µg/mL (100 µL/well) can bind Anti-Fusion glycoprotein F0 Antibody, Human IgG1 (D25) with a linear range of 0.2-1 ng/mL (QC tested).




Immobilized HRSV (A) Pre-fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H7) at 1 µg/mL (100 µL/well) on a Nickel Coated plate can bind Monoclonal Anti-HRSV F Antibody, Human IgG1 (AM14) with a linear range of 0.2-2 ng/mL (QC tested).




Immobilized HRSV (A) Pre-fusion glycoprotein F0, His Tag (MALS verified) (Cat. No. RSF-V52H7) at 2 µg/mL (100 µL/well) can bind Monoclonal Anti-RSV-F0 broadly Antibody, Mouse IgG1 (3B4) (Cat. No. RS0-Y161) with a linear range of 0.098-3.125 ng/mL (Routinely tested).



Immobilized HRSV (A) Pre-fusion glycoprotein F0, His Tag (MALS verified) (Cat. No. RSF-V52H7) at 2 µg/mL (100 µL/well) can bind Monoclonal Anti-RSV-F0 broadly Antibody, Mouse IgG1 (4G2) with a linear range of 0.024-0.781 ng/mL (Routinely tested).



Immobilized HRSV (A) Fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H7) at 2 µg/mL (100 µL/well) can bind Anti-Fusion glycoprotein F0 Antibody, Mouse IgG2a (101F) with a linear range of 0.2-5 ng/mL (Routinely tested).



Immobilized Anti-Fusion glycoprotein F0 Antibody, Mouse IgG2a (101F) at 2 µg/mL (100 µL/well) can bind HRSV (A) Fusion glycoprotein F0, His Tag (Cat. No. RSF-V52H7) with a linear range of 0.1-13 ng/mL (Routinely tested).

Discounts, Gifts,
and more!





Background

Human respiratory syncytial virus (HRSV) is the most common etiological agent of acute lower respiratory tract disease in infants and can cause repeated infections throughout life. The RSV fusion glycoprotein (RSV F) is the principal target of RSV neutralizing antibodies in human sera. The RSV F is a type I viral fusion protein synthesized as inactive, single-chain polypeptides that assemble into trimers. RSV F fuses the viral and host cell membranes by irreversible protein refolding from the labile prefusion conformation to the stable post-fusion conformation. Both states exhibit epitopes targeted by neutralizing antibodies, and post-fusion RSV F is being developed as a vaccine candidate.

Clinical and Translational Updates

Discounts, Gifts,
and more!

