

## Synonym

TACSTD2,GA733-1,M1S1,TROP2

#### **Source**

Cynomolgus TROP-2 Protein, His Tag(TR2-C52H5) is expressed from human 293 cells (HEK293). It contains AA Gln 31 - Thr 274 (Accession # XP 005543292.2).

Predicted N-terminus: Gln 31

#### **Molecular Characterization**

TROP-2(Gln 31 - Thr 274) XP\_005543292.2

Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 29.4 kDa. The protein migrates as 40-53 kDa when calibrated against <u>Star Ribbon Pre-stained Protein Marker</u> under reducing (R) condition (SDS-PAGE) due to glycosylation.

## Endotoxin

Less than 1.0 EU per µg by the LAL method.

## **Purity**

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

#### **Formulation**

Lyophilized from 0.22  $\mu m$  filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

#### Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

#### Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


## **SDS-PAGE**



Cynomolgus TROP-2 Protein, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With <u>Star Ribbon Pre-stained Protein Marker</u>).

# **Bioactivity-ELISA**

## **SEC-MALS**




The purity of Cynomolgus TROP-2 Protein, His Tag (Cat. No. TR2-C52H5) is more than 90% and the molecular weight of this protein is around 30-45 kDa verified by SEC-MALS.

Report

# Cynomolgus TROP-2 / TACSTD2 Protein, His Tag (MALS verified)







Mouse Monoclonal Antibody Against Human TROP-2 Protein, Mouse IgG1 Conc. (ng/mL)

Immobilized Cynomolgus TROP-2 Protein, His Tag (Cat. No. TR2-C52H5) at 1  $\mu$ g/mL (100  $\mu$ L/well) can bind Mouse Monoclonal Antibody Against Human TROP-2 Protein, Mouse IgG1 with a linear range of 0.1-2 ng/mL (QC tested).

## Background

TROP-2 is a single-copy gene in human cells, and encodes a type-1 transmembrane glycoprotein which is over-expressed in various malignancies, also referred to as tumor associated calcium signal transducer 2 (TACSTD2), GA733-1 or M1S1. TROP-2 is related to epithelial cell adhesion molecule (EpCAM), also called TROP-1, gp40, and KSA. Trop-1 and Trop-2 are homologous to serum IGF-II-binding proteins and appear as signal transducers. Thus, they likely represent novel cell-surface receptors and may play a role in regulating the growth of carcinoma cells.

# **Clinical and Translational Updates**

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

