

Synonym

TNFRSF10B, TRAILR2, TRAIL-R2, CD262, DR5, KILLER, TRICK2, ZTNFR9, TRICKB

Source

Cynomolgus TRAIL R2, His Tag (TR2-C52H7) is expressed from human 293 cells (HEK293). It contains AA Ile 58 - Ser 212 (Accession # [A0A2K5TXK0-1](#)). Predicted N-terminus: Ile 58

Molecular Characterization

DR5(Ile 58 - Ser 212) A0A2K5TXK0-1	Poly-his
---------------------------------------	----------

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 19.1 kDa. The protein migrates as 25-30 kDa under reducing (R) condition (SDS-PAGE) due to Glycosylation.

Endotoxin

Less than 1.0 EU per μg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4. Normally trehalose is added as protectant before lyophilization.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.


Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.


Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


SDS-PAGE

Cynomolgus TRAIL R2, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained overnight with Coomassie Blue. The purity of the protein is greater than 90%.

Bioactivity-ELISA**SEC-MALS**

The purity of Cynomolgus TRAIL R2, His Tag (Cat. No. TR2-C52H7) was more than 90% and the molecular weight of this protein is around 25-35 kDa verified by SEC-MALS.

[Report](#)



Immobilized Cynomolgus TRAIL R2, His Tag (Cat. No. [TR2-C52H7](#)) at 5 µg/mL (100 µL/well) can bind Human TRAIL, Mouse IgG2a Fc Tag (Cat. No. [TRL-H5259](#)) with a linear range of 2-63 ng/mL (QC tested).

Background

Tumor necrosis factor receptor superfamily member 10B (TNFRSF10B) is also known as TNF-related apoptosis-inducing ligand receptor 2 (TRAILR2), Death receptor 5 (DR5), CD262, KILLER, is a member of the TNF-receptor superfamily, and contains an intracellular death domain. TNFRSF10B / DR-5 is widely expressed in adult and fetal tissues; very highly expressed in tumor cell lines. TRAILR2 / CD262 / DR5 is the receptor for the cytotoxic ligand TNFSF10/TRAIL. The adapter molecule FADD (a death domain containing adaptor protein) of TRAIL-R2 / TNFRSF10B recruits caspase-8 to the activated receptor. The resulting death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis. CD262 / DR5 Promotes the activation of NF-kappa-B. DR5 is essential for ER stress-induced apoptosis and is regulated by p53/TP53.

References

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.